

Display Report


Analysis Info

Analysis Name D:\Data\test\0713\xinchengpiganerqingzhaertong_96_01_1916.d
Method 1916.m
Sample Name xinchengpiganerqingzhaertong
Comment

Operator
Instrument bruker
amaZon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	140.0 V	n/a	n/a	Trap Drive	50.6
Accumulation Time	2071 μ s	Averages	5 Spectra	Auto MS/MS	off



#	RT [min]	Area	Height	Area Frac. %
1	4.0	3.3448	0	0.1
2	5.6	2.6811	0	0.1
3	9.3	34.7101	2	0.8
4	11.1	3.4424	0	0.1
5	12.1	4176.6068	178	95.5
6	12.8	21.4345	1	0.5

Display Report

Analysis Info

Analysis Name D:\Data\test\0713\xinchengpiganerqingzhaertong_96_01_1916.d
Method 1916.m
Sample Name xinchengpiganerqingzhaertong
Comment

Operator bruce
Instrument amazon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	140.0 V	n/a	n/a	Map Drive	50.6
Accumulation Time	2071 μ s	Averages	5 Spectra	Auto MS/MS	off

