

Display Report

Analysis Info

Analysis Name D:\Data\TEST\1110\yiluoxinfugan.d
Method esi-1000.m
Sample Name yiluoxinfugan
Comment

Instrument amaZon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Negative	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	-140.0 V	n/a	n/a	Trap Drive	79.4
Accumulation Time	2215 μ s	Averages	5 Spectra	Auto MS/MS	off

